

Stacc-seq protocol (version 5)

Bofeng Liu

Xie Lab

March 23th, 2021

Stacc-seq: Small-scale Tn5-assisted Chromatin Cleavage with sequencing

Reagents

- Ethanol
- ddH₂O
- pA/G-Tn5 (Tn5 and proteinA/G fusion protein, Vazyme TD901, TD902)
- AMPure beads
- Concanavalin-coated magnetic beads (Polyscience, 86057)
- Qubit DNA HS
- Phenol-chloroform
- Carrier RNA (50x dilute first to make 20ng/ul solution)
- NaOAc (3M)
- Glycogen
- 5x TTBL (Vazyme, TD502)
- PCR kit for Tn5 (Vazyme, TD601)
- Digitonin (Thermo, BN2006)
- PCR primer for Tn5 (Vazyme, TD202)
- TE buffer (Tris-EDTA, for DNA)

- **Buffer1**

Tris-HCl (pH 7.4)	10mM
NaCl	150mM
Spermidine	0.5mM
Roche complete (EDTA free)	1X

- **Buffer2**

Tris-HCl (pH 7.4)	10mM
KCl	10mM
CaCl ₂	1mM
MnCl ₂	1mM

Materials

- Magnetic rack
- PCR thermal cycler
- Incubator
- Thermo mixer
- Low binding EP, PCR tube and pipette tube
- 4°C Centrifuge
- Phase-lock tube

Procedures

1. Harvest and count cells.

- Cells should be harvested freshly and frozen pellet cells do not guarantee good quality of Stacc-seq library
- Cells should be intact and in a homogenous, single-cell suspension.

2. Antibody and pA/G-Tn5 complex preparation

- Prepare D-buffer 1 (DB-1): add 1ul fully dissolved 5% digitonin (@95° C, >5min) into 1ml Buffer1 and vortex.
- Add 7ul DB-1, 0.5ul pA/G-Tn5, 0.5ug antibody into 200ul low-binding tube and vortex.
- Incubate at 4°C for 30min.

Stacc-seq without wash step:

3-1. Prepare cell lysis

- Add 1M~500 cells into 200ul low-binding tube (volumn < 1ul).
- Add 6ul DB-1 into cells and vortex.
- Incubate at 4°C for 10min, vortex each 2.5min for 3 times.

4-1. ChIC via antibody-Tn5 complex

- Add 35ul DB-1 into cell lysis, then add antibody and pA/G-Tn5 mixture.
- Add 12.5ul pre-warmed 5xTTBL.
- Incubate in Thermo mixer at 37°C for 30min (400rpm, 10s ON, 10s OFF)

Stacc-seq with wash step:

3-2. Prepare cell lysis

- Add 1M~500 cells into 200ul low-binding tube (volumn < 1ul).

- Suspend cell with 50ul Buffer 1.
- Add 10ul pre-washed ConA beads (wash with 200ul Buffer 2 once and resuspend in 10ul Buffer 2), and incubate at room temperature for 10min.
- Remove supernatant on a magnet stand, and then washed with 200ul Buffer 1 twice.
- Remove supernatant on a magnet stand, and then add pA/G–Tn5 mixture and DB-1 buffer (50 µl in total)
- Incubate on an Eppendorf Thermomixer at 4 °C, 400 rpm for 2 h.
- Remove supernatant on a magnet stand, and wash twice with 200ul DB-1 buffer.
- Suspend beads with 50ul DB-1 buffer.
- Add 12.5ul 5xTTBL.

4-2. ChIC via antibody-Tn5 complex

- Suspend beads with 50ul DB-1 buffer.
- Add 12.5ul pre-warmed 5xTTBL.
- Incubate in Thermo mixer at 37°C for 30min (400rpm, 10s ON, 10s OFF)

Follow step 4-1 or 4-2:

5. Phenol-chloroform Purification

- Add 2ul carrier RNA, 2ul spike-in DNA (optional) and 65ul TE into cell lysis.
- Transfer the mixture into a phase-lock tube.
- Add 130ul Phenol-chloroform and vortex for 15s.
- Incubate at room temperature for 5min.

- 12,000 rpm centrifuge 5 min and extract supernatant
- Add 650ul ethanol, 24ul NaOAc, 2ul glycogen
- Incubate at -20°C 30min or overnight (for low input cells).
- Spin at 4°C for 15 min at maximum speed.
- Remove supernatant and wash pellets with 75% ethanol. Air dry and resuspend pellets in 29ul H₂O.

6. PCR reaction

make the 50 µL PCR reaction mix in a 0.2 mL PCR tube to amplify transposed DNA fragments:

- Be careful to ensure that samples are barcoded appropriately for subsequent pooling and sequencing.
- The first 3min extension at 72 °C is critical to allow extension of both ends of the primer after transposition, thereby generating amplifiable fragments.

DNA product	29ul
5xTAB	10ul
N5XX	5ul
N7XX	5ul
TAE	1ul

Thermal cycle as follows:

72°C	3min
98°C	30sec
16cycle	98°C 15sec

	60°C	30sec
	72°C	3min
72°C		5min
4°C		hold

7. Purify amplified library using AMPure Beads.

- 0.4x beads were used to filter out large fragments
- 1.7x beads were used to filter out small fragments
- Elute the purified library in 20 μ L ddH₂O

8. Library QC and quantification

- Measure the concentration of libraries by using Qubit.
- Quantify relative amount of libraries by using qPCR Illumina protocol.
- Pool your libraries according to qPCR results above.